

EmosoftTM ointment

White Soft Paraffin 50% w/w
& Liquid Paraffin 50% w/w



Presentation

EmosoftTM ointment: Each gram ointment contains White Soft Paraffin BP 500 mg and Liquid Paraffin BP 500 mg.

Description

White Soft Paraffin and Liquid Paraffin are emollients with occlusive properties for topical application. There is no significant systemic absorption of the active ingredients but has substantial penetration into the stratum corneum.

Indications and Uses

For general use as an emollient in the symptomatic relief of dry skin conditions.

Dosage and Administration

Apply a thin film of the ointment to the affected area of the skin, in the direction of hair growth, as required. Repeat as necessary. Ideally the product should be applied three or four times a day or at least twice a day.

In adults, where a large area of the body is affected, up to 500 g a week may be used. For topical administration only.

Contraindications

Contraindicated to those patients with a history of sensitivity reactions to any of its components.

Side effects

Prolong use may cause folliculitis, should this occur, use of the product should be discontinued.

Precautions

The ointment should be applied in the direction of hair growth to reduce the incidence of folliculitis.

Use in pregnancy and lactation

The safety of White Soft Paraffin and Liquid Paraffin during pregnancy and lactation has not been established, but use during this periods is not considered to constitute a hazard.

Overdose

EmosoftTM ointment has low toxicity and if ingested treatment should be supportive and symptomatic. Gastric lavage should not be attempted following ingestion of paraffin-based products, due to the risk of inhalation or aspiration into the lungs.

Special precautions for storage

Do not store above 30⁰ C. Keep away from light and out of the reach of children.

Commercial pack

EmosoftTM Ointment: Each tube contains 50 gm ointment.



Manufactured by

Incepta Pharmaceuticals Ltd

Dhamrai Unit, Dhaka, Bangladesh

TM Trademark

EMT

V.N. 01